

Peptide mix

Mixture of two basics and one acidic peptide

Experimental

Column: SVEA C4, 5 μm 110 \AA 4.6 x 150 mm
Instrument: HPLC
Mobile phase: A: 0.1% TFA in water; B: 0.085% TFA in acetonitrile.
Gradient: 10-50% B in 60 minutes
Flow rate: 1.0 mL/min
Column temperature: 30°C
Detector: UV 210 nm
Sample: 1. 14 aminoacids pI 4.5
2. 9 aminoacids pI 11.7
3. 32 aminoacids pI 11.1

